



Recombinant *Yersinia pseudotuberculosis* serotype O:1b Maltoporin 1 (IamB1)

Product Code	CSB-EP412316YAN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7FDH5
Product Type	Recombinant Protein
Immunogen Species	<i>Yersinia pseudotuberculosis</i> serotype O:1b (strain IP 31758)
Purity	>85% (SDS-PAGE)
Sequence	VDFHGY ARSGIGWTAS GGEQQCFQTT GAQSKYRLGN ECETYAELKL GQELWKEGDK SFYLDTNVAY SVSQRDDWES TDPAFREANV QGKNLIESLP GSTIWAGKRF YQRHDVHMID FYYWDISGPG AGLETIDLGF GKLSVAATR SESGGSSAWI DNQRENAKYT INNVYDVRLA GLETNPGGSL ELGVDYGRAD TQEGYSLAPN ASKDGVMLTA EHTQSLMGGF NKFVVQYATD SMTSYNTGHS QGTSVNNNGH LLRVIDHGAI NLAEKWDMMY VALYQDIDLD NNNGNTWYSV GVRPMYKWTP IMSTLLEAGY DNVKSQHTGE RNGQYKLT LA QQWQAGDSIW SRPAIRVFAT YANWDEKWGY SDTTGVAQDG TIGTNSRGKN NEVTFGAQFE AWW
Source	E.coli
Target Names	IamB1
Protein Names	Recommended name: Maltoporin 1 Alternative name(s): Maltose-inducible porin 1
Expression Region	25-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.