



Recombinant Chicken F-actin-capping protein subunit beta isoforms 1 and 2 (CAPZB)

Product Code	CSB-BP004514CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P14315
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	SDQQLDCAL DLMRRLPPQQ IEKNLSDLID LVPSLCEDLL SSV DQPLKIA RDKVVGKDYL LCDYNRDGDS YRSPWSNKYD PPLEDGAMPS ARLRKLEVEA NNAFDQYRDL YFEGGVSSVY LWDL DHGFAG VILIKKAGDG SKKIKGCWDS IHVVEVQEKS SGRTAHYKLT STVMLWLQTN KTGSGTMNLG GSLTRQMEKD ETVSDSSPHI ANIGRLVEDM ENKIRSTLNE IYFGKTKDIV NGLRSIDAIP DNQKYKQLQR ELSQVLTQRQ IYIQPDN
Source	Baculovirus
Target Names	CAPZB
Protein Names	Recommended name: F-actin-capping protein subunit beta isoforms 1 and 2 Alternative name(s): Beta-actinin subunit II CapZ 36/32 CapZ B1 and B2
Expression Region	2-277
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	CAPZB is a member of the F-actin capping protein family. This gene encodes the beta subunit of the barbed-end actin binding protein. The protein regulates growth of the actin filament by capping the barbed end of growing actin filaments.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.