



Recombinant *Saccharomyces cerevisiae* Respiration factor 1 (RSF1)

Product Code	CSB-YP413820STA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6ZM92
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKDLNPEMGK FATTKGPPQD NRGMVDIATL PNFPANRSGT PREEMYLAPN KMETPRILNM NMVPDYLQKE NFSPDFSSAT VSAKSSPVNV THDESLPLGT IESNSTKDSK YAVQRQQQQV VDFIENNMQL LSSETLNFERS DIMKTLELPI PKRRDIKGNH LSKLLFAKSP LTINTYCQFY DRRTKRICNQ EMIWKDKNSR EKHGSRKYQR HLSKVHVDVQL TPNNFTEFFD HNSPLFQECY DYQSRLMRDL LVEPDAKFKE KKKKKKGDVN GNHPETGSSL INHQVQQQNV RELQSIAMN DLIEILIDLN IPFSVLDYQP MRNWLKYSI ISTDTPDEV YFKTDPGVNE LEHNSSNLNN SNSGTPHNNH QNQHTN
Source	Yeast
Target Names	RSF1
Protein Names	Recommended name: Respiration factor 1
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.