



# Recombinant *Saccharomyces cerevisiae* Adenylyltransferase and sulfurtransferase UBA4 (UBA4)

<b>Product Code</b>	CSB-MP413848STA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A6ZT19
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNDYHLEDTT SELEALRLEN AQLREQLAKR EDSSRDYPLS LEEYQRYGRQ MIVEETGGVA GQVKLKNTRV LVVGAGGLGC PALPYLAGAG VGQIGVDND VVETSNLHRQ VLHDSSRVGM LKCESARQYI TKLNPHINVV TYPVRLNSSN AFDIFKGYNY ILDCTDSPLT RYLVSDVAVN LGITVVSASG LGTEGQLTIL NFNNIGPCYR CFYPTPPPPN AVTSCQEGGV IGPCIGLVGT MMAVETLKLI LGIYTNENFS PFLMLYSGFP QQLRRTFKMR GRQEKCLCCG KNRTITKEAI EKGEINYELF CGARNYNVCE PDERISVDAF QRIYKDDEFL AKHIFLDVRP SHHYEISHFP EAVNIPIKNL RDMNGDLKKL QEKLPSVEKD SNIVILCRYG NDSQLATRLL KDKFGFSNVR DVRRGGYFKYI DDIDQTIPKY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	UBA4
<b>Protein Names</b>	Recommended name: Adenylyltransferase and sulfurtransferase UBA4 Alternative name(s): Needs CLA4 to survive protein 3 Ubiquitin-like protein activator 4 Including the following 2 domains: Adenylyltransferase UBA4 EC= 2.7.7.-
<b>Expression Region</b>	1-440
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.