



Recombinant *Saccharomyces cerevisiae* Spindle pole component 29 (SPC29)

Product Code	CSB-MP413864STA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6ZWC8
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MDYNNFGNSA SKKFQDDTLN RVRKEHEEAL KKLREENFSS NTSELGNKKH YRAQERMSSP LHRLSPAGKS DDRKVKSPLD DKLRRQLREG NTRLPPPPFS SYGMPPTNRS NLDRIRRRTS SPVRTDKFAS QNVIDDQRLE IKYLERIVYD QGTVIDNLTS RITRLESFIL NSISDRGDKN FASLEHSRSF SGFPTNKTYG LQMGGLYEND MPYRRSSDNI NKEGAREDRS SQIHIEENEST EDILKILSSS FHN
Source	Mammalian cell
Target Names	SPC29
Protein Names	Recommended name: Spindle pole component 29
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.