



Recombinant *Saccharomyces cerevisiae* tRNA (guanine-N (7)-)-methyltransferase

Product Code	CSB-EP413867STA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6ZXD2
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKAKPLSQDP GSKRYAYRIN KEENRKELKH VKINESSLVQ EGQKIDLPKK RYYRQRAHSN PFSHQLEYP VSPQDMDWSK LYPYYKNAEN GQMTKKVTIA DIGCGFGGLM IDLSPAPPED LILGMEIRVQ VTNYVEDRII ALRNNTASKH GFQINIVLRG NAMKFLPNFF EKGQLSKMFF CFPDPHFKQR KHKARIITNT LLSEYAYVLK EGGVVYTITD VKDLHEWMVK HLEEHPLEFER LSKEWEENDE CVQIMRNATE EGKKVERKKG DKFVACFTRL PTPAIL
Source	<i>E.coli</i>
Target Names	TRM8
Protein Names	Recommended name: tRNA (guanine-N(7)-)-methyltransferase EC= 2.1.1.33 Alternative name(s): Transfer RNA methyltransferase 8 tRNA (guanine(46)-N(7))-methyltransferase tRNA(m7G46)-methyltransferase
Expression Region	1-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.