



Recombinant *Yersinia pseudotuberculosis* serotype O:1b Maltoporin 2 (IamB2)

Product Code	CSB-YP413937YAN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7FF76
Product Type	Recombinant Protein
Immunogen Species	<i>Yersinia pseudotuberculosis</i> serotype O:1b (strain IP 31758)
Purity	>85% (SDS-PAGE)
Sequence	IEKIDFH GYMRAVGVS SDGGLAEWQK TMVGR LGNES DTYGEIGLGA EVYKKEDVSF YLESMVSM LS DGSNDSETTI GDDAQFGLRQ LNLQIKGLIP GDKEAVIWGG KRYRQRHDLH IIDTKYWNIS GSGAGIENYT VGPGAVSVAW VRGDANDVDT RITGDSDVNI NYIDVRYAGF KPWAGSWTEV GIDYAMPNPT KQQKEYGGLY DADNAVMLTG EISQDMFGGY NKLVLQYANK GLAQN MISQG GGWYDMWHKT DEAKGYRVIN TGLIPITDKF SFNHVLTWGS ANDITEYTDK TNLISLVGRA QYQFTQYVRA IGEVGGFYQK DTYHNGSNYK QGGEKYTIAL GLAEGPDFLS RPELRVFASY LNDS ENKPF EDGTSNDTWN FGVQVEAWW
Source	Yeast
Target Names	IamB2
Protein Names	Recommended name: Maltoporin 2 Alternative name(s): Maltose-inducible porin 2
Expression Region	24-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.