



Recombinant Mouse Granzyme A (Gzma)

Product Code	CSB-BP010081MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11032
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	II GGDTVPHSR PYMALLKLSS NTICAGALIE KNWVLTAHC NVGKRSKFIL GAHSINKEPE QQILTVKKAF PYPCYDEYTR EGDQLVRLK KKATVNRNVA ILHLPKKGDD VKPGTRCRVA GWGRFGNKSA PSETLREVNI TVIDRKICND EKHYNFHPVI GLNMICAGDL RGGKDSCNGD SGSPLLCDGI LRGITSFGGE KCGDRRWPGV YTFLLSDKHLN WIKKIMKGSV
Source	Baculovirus
Target Names	Gzma
Protein Names	Recommended name: Granzyme A EC= 3.4.21.78 Alternative name(s): Autocrine thymic lymphoma granzyme-like serine protease CTLA-3 Fragmentin-1 T cell-specific serine protease 1 Short name= TSP-1
Expression Region	29-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Cytolytic T lymphocytes (CTL) and natural killer (NK) cells share the remarkable ability to recognize, bind, and lyse specific target cells. They are thought to protect their host by lysing cells bearing on their surface nonself antigens, usually peptides or proteins resulting from infection by intracellular pathogens. The protein described here is a T cell- and natural killer cell-specific serine protease that may function as a common component necessary for lysis of target cells by cytotoxic T lymphocytes and natural killer cells.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.