



Recombinant Bovine Galectin-1 (LGALS1)

Product Code	CSB-EP012882BO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11116
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	ACGLVASNL NLKPGECLRV RGEVAADAKS FLLNLGKDDN NLCLHFNPRF NAHGVDVNTIV CNSKDAGAWG AEQRESAFPFP QPGSVVEVCI SFNQDTLTIK LPDGYEFKFP NRLNLEAINY LSAGGDFKIK CVAFE
Source	E.coli
Target Names	LGALS1
Protein Names	Recommended name: Galectin-1 Short name= Gal-1 Alternative name(s): 14 kDa lectin Beta-galactoside-binding lectin L-14-I Galaptin Lactose-binding lectin 1 Lectin galactoside-binding soluble 1 S-Lac lectin 1
Expression Region	2-135
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. This gene product may act as an autocrine negative growth factor that regulates cell proliferation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.