



Recombinant Sheep Deoxyribonuclease-1 (DNASE1)

Product Code	CSB-EP007049SH-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11937
Product Type	Recombinant Protein
Immunogen Species	Ovis aries (Sheep)
Purity	≥85% (SDS-PAGE)
Sequence	LKIAAFNIRT FGETKMSNAT LSSYIVRILR RYDIALIEQV RDSHLVAVGK LLDDLNQDDP NSYHYVVSEF LGRNSYKERY LFFRPNKVS VLDTYQYDDG CESCGNDSFS REPAVVKFSS PSTKVKAFI VPLHSAPSDA VAEINSLYDV YLDVQQKWDL NDIMLMGDFN ADCSYVTSSQ WSSIRLRTSS TFQWLIPDSA DTTATSTNCA YDRIVVAGSL LQSSVVGPSA VPFDFQAAYG LSNEMALAI DHYPVEVTLT
Source	E.coli
Target Names	DNASE1
Protein Names	Recommended name: Deoxyribonuclease-1 EC= 3.1.21.1 Alternative name(s): Deoxyribonuclease I Short name= DNase I
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes a member of the DNase family. This protein is stored in the zymogen granules of the nuclear envelope and functions by cleaving DNA in an endonucleolytic manner. At least six autosomal codominant alleles have been characterized, DNASE1*1 through DNASE1*6, and the sequence of DNASE1*2 represented in this record. Mutations in this gene have been associated with systemic lupus erythematosus (SLE), an autoimmune disease. A recombinant form of this protein is used to treat the one of the symptoms of cystic fibrosis by hydrolyzing the extracellular DNA in sputum and reducing its viscosity. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.