



Recombinant Vanderwaltozyma polyspora Autophagy-related protein 3 (ATG3)

| | |
|--------------------------|---|
| Product Code | CSB-BP419267VDW |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | A7TK16 |
| Product Type | Recombinant Protein |
| Immunogen Species | Vanderwaltozyma polyspora (strain ATCC 22028 / DSM 70294) (Kluyveromyces polysporus) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MLRSTLSSWR EYLTPVSHKS TFLSSGQITP DEFVQAGDYL CHMFPTWEWN KAGNDVLFNRN FLPEDKQFLV MRKVPCNVRA KQFVNIDDSA SEAFVQGIND EDGSIENGWM IGGGETDHL S KHSLSGSDVT PSGNNTKTMD ETPDIDDLM ESMGIEDTDD IIDTPQNTDR RYYDLYITYS TSYKVPKMYI VGFNGSGSPL TPEEMFEDIA PDYRSKTATI EKLPFYKNTI SSVSIHPCKH ANVMRILLDK VLVVKRRRRE EMSENHNEHK PNPESEDEWED LQNDVDDTIR ADQYLIIFLK FITSVTPGIE HDYTMEGW |
| Source | Baculovirus |
| Target Names | ATG3 |
| Protein Names | Recommended name: Autophagy-related protein 3 Alternative name(s): Autophagy-related E2-like conjugation enzyme ATG3 |
| Expression Region | 1-318 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |