



Recombinant Escherichia coli O139:H28 Autoinducer 2-binding protein IsrB (IsrB)

Product Code	CSB-MP419491EJD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZLX3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Purity	≥85% (SDS-PAGE)
Sequence	AERI AFIPKLVGVG FFTSGGNGAQ QAGKELSVDV TYDGPTEPSV SGQVQLINNF VNQGYNAIIV SAVSPDGLCP ALKRAMQRGV RVLTWSDTK PECRSYYINQ GTPAQLGGML VDMAARQVNK DKAKVAFFYS SPTVTDQNQW VKEAKAKIAK EHPGWEIVTT QFGYNDATKS LQTAEGILKA YSDLDAIAP DANALPAAAQ AAENLKNKDV AIVGFSTPNV MRPYVERGTV KEFGLWDVVQ QGKISVYVAD ALLKKGSMKT GDKLDIQGVG QVEVSPNSVQ GYDYEADGNG IVLLPERVIF NKENIGKYDF
Source	Mammalian cell
Target Names	IsrB
Protein Names	Recommended name: Autoinducer 2-binding protein IsrB Short name= AI-2-binding protein IsrB
Expression Region	27-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.