



Recombinant Escherichia coli O139:H28 3-phenylpropionate/cinnamic acid dioxygenase ferredoxin subunit (hcaC)

Product Code	CSB-YP419512EJD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZPY3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Purity	>85% (SDS-PAGE)
Sequence	MNRIYACPVA DVPEGEALRI DTSPVIALFN VGGEFYAIND RCSHGNASMS EGYLEDDATV ECPLHAASFC LKTGKALCLP ATDPLTTYPV HVEGGDIFID LPEAQP
Source	Yeast
Target Names	hcaC
Protein Names	Recommended name: 3-phenylpropionate/cinnamic acid dioxygenase ferredoxin subunit Alternative name(s): Digoxigenin ferredoxin subunit
Expression Region	1-106
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.