



Recombinant Escherichia coli O139:H28 ATP synthase gamma chain (atpG)

Product Code	CSB-MP419540EJD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZTU5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Purity	>85% (SDS-PAGE)
Sequence	MAGAKEIRSK IASVQNTQKI TKAMEMVAAS KMRKSQDRMA ASRPYAETMR KVIGHLAHGN LEYKHPYLEE RDVKRVGYLV VSTDRLCGG LNINLFKLL AEMKTWTDKG VQCDLAMIGS KGVSVFFNSVG GNVVAQVTGM GDNPSLSELI GPVKVMLQAY DEGRDKLYI VSNKFINTMS QVPTISQLLP LPASDDDDLK HKSWDYLYEP DPKALLDTLL RRYVESQVYQ GVVENLASEQ AARMVAMKAA TDNGGSLIKE LQLVYNKARQ ASITQELTEI VSGAAAV
Source	Mammalian cell
Target Names	atpG
Protein Names	Recommended name: ATP synthase gamma chain Alternative name(s): ATP synthase F1 sector gamma subunit F-ATPase gamma subunit
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.