



Recombinant Escherichia coli O9:H4 Bifunctional protein FoID (foID)

Product Code	CSB-YP419569EJF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZXI3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	>85% (SDS-PAGE)
Sequence	MAAKIIDGKT IAQQVRSEVA QKVQARIAAG LRAPGLAVVL VGSNPASQIY VASKRKACEE VGFVRSYDL PETTSEAELL ELIDALNADN TIDGILVQLP LPAGIDNVKV LERIHDPKDV DGFHPYVGR LCQRAPRLRP CTPRGIVTLL ERYNIDTFGL NAVVIGASNI VGRPMSMELL LAGCTTTVTH RFTKNLRHHV ENADLLIVAV GKPGFIPGDW IKEGAIVIDV GINRLENGKV VGDVVFEDAA KRASYITPVP GGVGPMTVAT LIENTLQACV EYHDPQGE
Source	Yeast
Target Names	foID
Protein Names	Recommended name: Bifunctional protein FoID Including the following 2 domains: Methylenetetrahydrofolate dehydrogenase EC= 1.5.1.5 Methenyltetrahydrofolate cyclohydrolase EC= 3.5.4.9
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.