



# Recombinant Escherichia coli O139:H28 Pantothenate synthetase (panC)

<b>Product Code</b>	CSB-EP419882EJD-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A7ZHM3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O139:H28 (strain E24377A / ETEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLIETLPLL RQQIRRLRME GKRVALVPTM GNLHDGHMKL VDEAKARADV VVVSIFVNPM QFDRPEDLAR YPRTLQEDCE KLNKRKVDLV FAPSVKEIYP NGTETHTYVD VPGLSTMLEG ASRPGHFRGV STIVSKLFNL VQPDIACFGE KDFQQLALIR KVMADMGFDI EIVGVPIMRA KDGLALSSRN GYLTAEQRKI APGLYKVLSS IADKLQAGER DLDEIITIAG QELNEKGFRA DDIQIRDADT LLEVSETSKR AVILVAAWLG DARLIDNKMV ELA
<b>Source</b>	E.coli
<b>Target Names</b>	panC
<b>Protein Names</b>	Recommended name: Pantothenate synthetase Short name= PS EC= 6.3.2.1 Alternative name(s): Pantoate--beta-alanine ligase Pantoate-activating enzyme
<b>Expression Region</b>	1-283
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.