



Recombinant Escherichia coli O139:H28 Hydroxylamine reductase (hcp)

Product Code	CSB-BP419898EJD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZJU2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Purity	>85% (SDS-PAGE)
Sequence	MFCVQCEQTI RTPAGNGCSY AQGMCGKTAE TSDLQDLLIA ALQGLSAWAV KAREYGIINH DVDSFAPRAF FSTLTNVNFD SPRIVGYARE AIALREALKA QCLAVDANAR VDNPMANLQL VSDDLGE LQR QAAEFTP NKD KAAIGENILG LRLLC LYG LK GAAAYMEHAH VLGQYDNDIY AQYHKIMAWL GTWPADMNAL LECSMEIGQM NFKVMSILDA GETGKYGHPT PTQVNVKATA GKCILISGHD LKDLYNLLEQ TEGTGVNVYT HGEMLP AHGY PELRKFKHLV GNYGSGWQNQ QVEFARFPGP IVMTSNCIID PTVGAYDDRI WTRSIVGWPG VRHLDGEDFS AVIAQAQQMA GFPYSEIPHL ITVGFG RQTL LGAADTLIDL VSREKLRHIF LLGGCDGARG ERHYFTDFAT SVPDDCLILT LACGKYRFNK LEFGDIEGLP RLVDAGQCND AYSAILAVT LAEKLGC GVN DLPLSLVLSW FEQKAIVILL TLLSLGVKNI VTGPTAPGFL TPDLLAVLNE KFGLRSITTV EEDMKQLLSA
Source	Baculovirus
Target Names	hcp
Protein Names	Recommended name: Hydroxylamine reductase EC= 1.7.-.- Alternative name(s): Hybrid-cluster protein Short name= HCP
Expression Region	1-550
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.