



# Recombinant Escherichia coli O139:H28 Periplasmic trehalase (treA)

<b>Product Code</b>	CSB-MP419906EJD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A7ZKW9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O139:H28 (strain E24377A / ETEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	EETPVTPQPP DILLGPLFND VQNAKLFPDQ KTFADAVPNS DPLMILADYR MQQNQSGFDL RHFVNVNFTL PKEGEKYVPP EGQSLREHID GLWPVLTRST ENTEKWDSLL PLPEPYVVP GGRFREYVYWD SYFTMLGLAE SGHWDKVADM VANFAHEIDT YGHIPNGNRS YYLSRSQPPF FAMMVELLAQ HEGDAALKQY LPQMKEYAY WMDGVENLQA GQQEKRIVKL QDGTLLNRYW DDRDTPRPES WVEDIATAKS NPNRPATEIY RDLRSAASG WDFSSRWMDN PQQNLNLTTRT SIVPVDLNSL MFKMEKILAR ASKAAGDNAM ANQYETLANA RQKGIEKYLW NDQQGWYADY DLKSHKVRNQ LTAAALFPLY VNAAAKDRAN KMATATKTHL LQPGGLNTTS VKSGQQWDAP NGWAPLQWVA TEGLQNYGQK EVAMDISWHF LTNVQHTYDR EKKLVEKYDV STTGTGGGGG EYPLQDGF GW TNGVTLKMLD LICPKEQPCD NVPATRPVK SATTQPSTKE AQPTP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	treA
<b>Protein Names</b>	Recommended name: Periplasmic trehalase EC= 3.2.1.28 Alternative name(s): Alpha,alpha-trehalase Alpha,alpha-trehalose glucohydrolase
<b>Expression Region</b>	31-565
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.