



Recombinant Escherichia coli O9:H4 Peptide chain release factor 2 (prfB)

Product Code	CSB-EP420026EJF-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8A432
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	>85% (SDS-PAGE)
Sequence	MFEINPVNNR IQDLTERS DV LRGYLCYDAK KERLEEVNAE LEQPDVWNEP ERAQALGKER SSLEAVVDTL DQMKQGLEDV SGLLELAVEA DDEETFNEAV AELDALEEKL AQLEFRRMFS GEYDSADCYL DIQAGSGGTE AQDWASMLER MYLRWAESRG FKTEIIIESE GEVAGIKSVT IKISGDYAYG WLRTETGVHR LVRKSPFDSG GRRHTSFSSA FVYPEVDDDI DIEINPADLR IDVYRASGAG GQHVNRTESA VRITHIPTGI VTQCQNDRSQ HKNKDQAMKQ MKAKLYELEM QKKNAEKQAM EDNKSDIGWG SQIRSYVLDD SRIKDLRTGV ETRNTQAVLD GSLDQFIEAS LKAGL
Source	E.coli
Target Names	prfB
Protein Names	Recommended name: Peptide chain release factor 2 Short name= RF-2
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.