



Recombinant Escherichia coli O139:H28 UDP-2,3-diacylglucosamine hydrolase (IpxH)

Product Code	CSB-YP420905EJD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZIT3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Purity	>85% (SDS-PAGE)
Sequence	MATLFIADLH LCVEEPAITA GFLRFLAGEA RKADALYILG DLFEAWIGDD DPNPLHRQMA AAIKAVSDSG VPCYFIHG NR DFLLGKRFAR ESGMTLLPEE KVLELYGRRV LIMHGDTLCT DDAGYQAFRA KVHKPWLQML FLALPLFVRK RIAARMRANS KEANSSKSLA IMDVNQNAV V SAMEKHQVQW LIHGHTHRPA VHELIANQQP AFRVVLGAWH TEGSMVKVTA DDVELIHFPF
Source	Yeast
Target Names	IpxH
Protein Names	Recommended name: UDP-2,3-diacylglucosamine hydrolase EC= 3.6.1.54 Alternative name(s): UDP-2,3-diacylglucosamine diphosphatase
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.