



Recombinant Escherichia coli O139:H28 Soluble pyridine nucleotide transhydrogenase (sthA)

Product Code	CSB-MP420978EJD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZUI2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Purity	>85% (SDS-PAGE)
Sequence	<p>MPHSYDYDAIVIGSGPGGEGAAMGLVKQGARVAVIERYQNVGGGCTHWGTIP SKALRHAV SRIIEFNQNPLYSDHSRLLRSSFADILNHADNVINQQTRMRQGFYERNHCEILQ GNARFV DEHTLALDCPDGSVETLTAEKFVIACGSRPYHPTDVFTHPRIYDSDSILSMHH EPRHVL IYGAGVIGCEYASIFRGMVDKVDLINTRDLLAFLDQEMSDSLSYHFWNSGVVI RHNEEY EKIEGCDDGVIMHLKSGKCLKADCLLYANGRTGNTDSLALQNIGLETDSRGQL KVNSMYQ TAQPHVYAVGDVIGYPSLASAAYDQGRIAAQALVKGEATAHLIEDIPTGIYTIPEI SSVG KTEQQLTAMKVPYEVGRAQFKHLARAQIVGMNVGTLKILFHRETKEILGIHCFG ERAAEI IHIGQAIMEQKGGGNTIEYFVNTTFNYPTMAEAYRVAALNGLNRLF</p>
Source	Mammalian cell
Target Names	sthA
Protein Names	Recommended name: Soluble pyridine nucleotide transhydrogenase Short name= STH EC= 1.6.1.1 Alternative name(s): NAD(P)(+) transhydrogenase [B-specific]
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.