



Recombinant Escherichia coli O9:H4 Alkanesulfonate monooxygenase (ssuD)

Product Code	CSB-EP421011EJF-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZYN8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	≥85% (SDS-PAGE)
Sequence	MSLNMFWFLP THGDGHYLG TEEGSRPVDHG YLQQIAQAAD RLG YTGVLIP TGRSCEDAWL VAASMIPVTQ RLKFLVALRP SVTSPTVAAR QAATLDRLSN GRALFNLVTG SDPQELAGDG VFLDHSEYEA ASAEFTQVWR RLLQRETVDF NGKHIHVRGA KLLFPAIQQP YPPLYFGGSS DVAQELAAEQ VDL YLTWGEP PELVKEKIEQ VRAKAAAHGR KIRFGIRLHV IVRETNDEAW QAAERLISHL DDETIKAQA AFARTDSVGQ QRMAALHNGK RDNLEISPNL WAGVGLVRGG AGTALVGDGP TVAARINEYA ALGIDSFVLS GYPHLEEAYR VGELLFPLLD VAIPEIPQPQ PLNPQGEAVA NDFIPRKVAQ S
Source	E.coli
Target Names	ssuD
Protein Names	Recommended name: Alkanesulfonate monooxygenase EC= 1.14.14.5 Alternative name(s): FMNH ₂ -dependent aliphatic sulfonate monooxygenase
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.