



# Recombinant Escherichia coli O9:H4 Autoinducer 2 import ATP-binding protein LsrA (IsrA)

<b>Product Code</b>	CSB-EP421021EJF-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A8A066
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O9:H4 (strain HS)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQTS DTRALP LLCARSVYKQ YSGVNV LKGI DFTLHQGEVH ALLGGNGAGK STLMKIIAGI TPADSGTLEI GGNNYARLTP VHAHQ LGIYL VPQEPLLFSS LSIKENILFG LAKKQLSMQK MKNLLAALGC QFDLHSLAGS LDVADRQMVE ILRGLMRDSR ILILDEPTAS LTPAETERLF TRLQELLATG VGIVFISHKL PEIRQIADRI SVMRDGTIAL SGKTSELSTD DIIQAITPVV REKSL SASQK LWLELPGNRP QHAAGTPVLT LENLTGEGFR NVSLTLNAGE ILGLAGLVGA GRTELAETLY GLRTL RGGRI MLNGNEINKL STGERLLRGL VYLPEDRQSS GLNLDASLAW NVCALTHNLR GFWAKTAKDN ATLERYRRAL NIKFNQPEQA ARTLSGGNQQ KILIAKCLEA SPQVLIVDEP TRGVDVSARN DIYQLLR SIA AQNVAVLLIS SDLEEIELMA DRVYVMHQGE IAHSALTGRD INVETIMRVA FGDSQRQEAS C
<b>Source</b>	E.coli
<b>Target Names</b>	IsrA
<b>Protein Names</b>	Recommended name: Autoinducer 2 import ATP-binding protein LsrA Short name= AI-2 import ATP-binding protein LsrA EC= 3.6.3.-
<b>Expression Region</b>	1-511
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.