



Recombinant Escherichia coli O9:H4 Arginine N-succinyltransferase (astA)

Product Code	CSB-YP421028EJF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8A0T9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	>85% (SDS-PAGE)
Sequence	MMVIRPVERS DVSALMQLAS KTGGLTSLP ANEATLSARI ERAIKTWQGE LPKSEQGYVF VLEDSETGTV AGICAIEVAV GLNDPWYNYR VGTLVHASKE LNVYNALPTL FLSNDHTGSS ELCTLFLDPD WRKEGNGYLL SKSRFMFMAA FRDKFNDKVV AEMRGVIDEH GYSPFWQSLG KRFFSMDFSR ADFLCGTGQK AFIAELMPKH PIYTHFLSQE AQDVIQVHP QTAPARAVLE KEGFRYRNYI DIFDGGPTLE CDIDRVRAIR KSRLVEVAEG QPAQGDFFAC LVANENYHHF RVVLRVTDPA TERLILTAAQ LDALKCHAGD RVRLVRLCAE EKTA
Source	Yeast
Target Names	astA
Protein Names	Recommended name: Arginine N-succinyltransferase Short name= AST EC= 2.3.1.109 Alternative name(s): AOST
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.