



# Recombinant Escherichia coli O9:H4 Bifunctional polymyxin resistance protein ArnA (arnA)

<b>Product Code</b>	CSB-BP421037EJF
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A8A2C2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O9:H4 (strain HS)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKTVVVFAYHD MGCLGIEALL AAGYEISAIF THTDNPGEKA FYGSVAHLAA ERGIPVYAPD NVNHPLWVER IAQLSPEVIF SFYYRHLICD EILQLAPAGA FNLHGSELLPK YRGRAPLNWV LVNGETETGV TLHRMVKRAD AGAIVAQLRV AIAPDDIAIT LHHKLCHAAR QLLEQTLPAI KHGNILEIAQ RENEATCFGR RTPDDSFLEW HKPASVLHNM VRAVADPWPG AFSYVGNQKF TVWSSRVHPH ASKAQPGSVI SIAPLLIACG DGALEIVTGQ AGDGITMQGS QLAQTGLLVQ GSRLNSQPAC TARRRTRVLI LGVNGFIGNH LTERLLREDH YEVYGLDIGS DAISRFLNHP HFHFVEGDIS IHSEWIEYHV KKCDVVLPLV AIATPIEYTR NPLRVFELDF EENLRIIRYC VKYRKRIIFP STSEVYGMCS DKYFDEDHSN LIVGPV NKPR WIYSVSKQLL DRVIWAYGEK EGLQFTLFRP FNWMGPRLDN LNAARIGSSR AITQLILNLV EGSPIKLIDG GKQKRCFTDI RDGIEALYRI IENAGNRCDG EIINIGNPEN EASIEELGEM LLASFEKHPL RHHFPPFAGF RVVESSSYG KGYQDVEHRK PSIRNAHRCL DWEPKIDMQE TIDETLDFFL RTVDLTDKPS
<b>Source</b>	Baculovirus
<b>Target Names</b>	arnA
<b>Protein Names</b>	Recommended name: Bifunctional polymyxin resistance protein ArnA Including the following 2 domains: UDP-4-amino-4-deoxy-L-arabinose formyltransferase EC= 2.1.2.13 Alternative name(s): ArnAFT UDP-L-Ara4N formyltransferase UDP
<b>Expression Region</b>	1-660
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.