



# Recombinant Escherichia coli O9:H4 L-fucose isomerase (fucl)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP421051EJF   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | A8A3T7  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O9:H4 (strain HS)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MKKISLPKIG IRPVIDGRRM GVRESLEEQT MNMAKATAAL LTEKLRHACG<br>AAVECVISDT CIAGMAEAAA CEEKFSSQNV GLTITVTPCW CYGSETIDMD<br>PTRPKAIWGF NGTERPGAVY LAAALAAHSQ KGIPAFSIYG HDVQDADDT<br>IPADVEEKLL RFARAGLAVA SMKGKSYLSL GGVSMGIAGS IVDHNFESW<br>LGMKVQAVDM TELRRRIDQQ IYDEAELEMA LAWADKNFRY GEDENNKQYQ<br>RNAEQSRAVL RESLLMAMCI RDMMQGNL ADIGRVEESL GYNAIAAGFQ<br>GQRHWTDQYP NGDTAEAILN SSFDWNGVRE PFVVATENDS<br>LNGVAMLMGH QLTGTAQVFA DVRTYWSPEA IERVTGHKLD GLAEHGIHL<br>INSGSAALDG SCKQRDSEGN PTMKPHWEIS QQEADACLAA TEWCPAIEY<br>FRGGGYSSRF LTEGGVPFTM TRVNIKGLG PVLQIAEGWS VELPKDVHDI<br>LNKRTNSTWP TTWFAPRLTG KGPFTDVYSV MANWGANHGV LTIGHVGADF<br>ITLASMLRIP VCMHNVEETK VYRPSAWAAH GMDIEGQDYR ACQNYGPLYK<br>R |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | fucl  |
| <b>Protein Names</b>     | Recommended name: L-fucose isomerase EC= 5.3.1.25 Alternative name(s): 6-deoxy-L-galactose isomerase Fuclase  |
| <b>Expression Region</b> | 1-591   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,   |



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.