



# Recombinant Human Ski oncogene (SKI)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP021356HU  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P12755  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MEAAAGGRGC FQPHPGLQKT LEQFHLSSMS SLGGPAAFSA<br>RWAQEAYKKE SAKEAGAAAV PAPVPAATEP PPVHLPAIQ PPPPVLP<br>FMPDRSTER CETVLEGETI SCFVVGGEKR LCLPQILNSV LRDFSLQ<br>AVCDELHIYC SRCTADQLEI LKVMGILPFS APSCGLITKT DAERLCN<br>YGGAYPPPCK KELAASLALG LELSERSVRV YHECFGKCKG LLVPELY<br>SAACIQCLDC RLMYPPHKFV VSHKALENR TCHWGFDSAN WRAYILLSQ<br>YTGKEEQARL GRCLDDVKEK FDYGNKYKRR VPRVSSEPPA SIRPKTD<br>SQSPAPSEKD KPSSWLRTLA GSSNKS LGCV HPRQRLSAFR PWSPAVS<br>KELSPHLPAL IRDSFY SYKS FETAVAPNVA LAPPAAQKV VSSPPCAA<br>RAPEPLATCT QPRKRKLTVD TPGAPETLAP VAAPEEDKDS EAEVEVE<br>EFTSSLSSLS SPSFTSSSSA KDLGSPGARA LPSAVPDAAA PADAPSG<br>ELEHLRQALE GGLDTKEAKE KFLHEVVKMR VKQEEKLSAA LQAKRSL<br>LEFLRVAKKE KLREATEAKR NLRKEIERLR AENEKMKKEA NESRLRL<br>LEQARQARVC DKGCEAGRRLR AKYSAQIEDL QVKLQHAHAD REQLR<br>EREAREHLEK VVKELQEQLW PRARPEAAGS EGAAELEP |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | SKI   |
| <b>Protein Names</b>     | Recommended name: Ski oncogene Alternative name(s): Proto-oncogene c-Ski  |
| <b>Expression Region</b> | 1-728   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full length protein   |
| <b>Target Details</b>    | This gene encodes the nuclear protooncogene protein homolog of avian sarcoma viral (v-ski) oncogene. It functions as a repressor of TGF-beta signaling, and may play a role in neural tube development and muscle differentiation.  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,   |



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.