



Recombinant Escherichia coli O139:H28 Betaine aldehyde dehydrogenase (betB)

Product Code	CSB-YP421920EJD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZI51
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Purity	>85% (SDS-PAGE)
Sequence	MSRMAEQQLY IHGGYTSATS GRTFETINPA NGNVLATVQA AGREDVDRAV KSAQQGQKIW AAMTAMERSR ILRRAVDILR ERNDELAKLE TLDTGKAYSE TSTVDIVTGA DVLEYYAGLI PSLEGSQIPL RETSFVYTRR EPLGVVAGIG AWNYPPIQIAL WKSAPALAAG NAMIFKPSEV TPLTALKLAE IYSEAGLPDG VFNVLPGVGA ETGQYLTEHP GIAKVSFTGG VASGKKVMAN SAASSLKEVT MELGGKSPLI VFDDADLDLA ADIAMMANFF SSGQVCTNGT RVFVPAKCKA AFEQKILARV ERIRAGDVFD PQTNFGPLVS FPHRDNVRLRY IVKGKEEGAR VLCGGDVLKG DGLDNGAWVA PTVFTDCSDE MTIVREEIFG PVMSILTYES EEEVIRRAND TDYGLAAGIV TADLNRAHRV IHQLEAGICW INTWGESPAE MPVGGYKHSI IGRENGVMTL QSYTQVKSII VEMAKFQSIF
Source	Yeast
Target Names	betB
Protein Names	Recommended name: Betaine aldehyde dehydrogenase Short name= BADH EC= 1.2.1.8
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.