



# Recombinant Escherichia coli O139:H28 Phosphate acyltransferase (plsX)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP421937EJD  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | A7ZKJ3   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli O139:H28 (strain E24377A / ETEC)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MTRLTLALDV MGGDFGPSVT VPAALQALNS NSQLTLLLVG NPDAITPLLA<br>KADFEQRSRL QIIPAQSVIA SDARPSQAIR ASRGSSMRVA LELVKEGRAQ<br>ACVSAGNTGA LMGLAKLLLK PLEGIERPAL VTVLPHQQKG KTVVLDLGAN<br>VDCDSTMLVQ FAIMGSVLAE EVVEIPNPRV ALLNIGEEEV KGLDSIRDAS<br>AVLKTIPSIN YIGYLEANEL LTGKTDVLC DGFTGNVTLK TMEGVVRMFL<br>SLLKSQGEGK KRSWWLLLLK RWLQKSLTRR FSHLNPQYN GACLLGLRGT<br>VIKSHGAANQ RAFAVAIEQA VQAVQRQVPQ RIAARLESVY PAGFELLDGG<br>KSGTLR |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | plsX   |
| <b>Protein Names</b>     | Recommended name: Phosphate acyltransferase EC= 2.3.1.n2 Alternative name(s): Acyl-ACP phosphotransacylase Acyl-[acyl-carrier-protein]--phosphate acyltransferase Phosphate-acyl-ACP acyltransferase   |
| <b>Expression Region</b> | 1-356  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |