



# Recombinant Escherichia coli O139:H28 Nicotinate-nucleotide--dimethylbenzimidazole phosphoribosyltransferase (cobT)

<b>Product Code</b>	CSB-BP421952EJD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A7ZNF5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O139:H28 (strain E24377A / ETEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQTLADLLNT IP AIDPAAMS RAQRHIDGLL KPVGSLGRLE ALAIQLAGMP GLNGIPHV GK KAVLVMCADH GVWEEGVAIS PKEVTAIQAE NMTRGTTGVC VLAAQAGANV YVIDVGIDTA EPIPGLINMR VARGSGNIAS APAMSRHQAE KLLLDVICYT RELAKNGVTL FVGELGMAN TTPAAAI VST ITGRAPEEVV GIGANLPTDK LANKIDV VRR AITLNQPNPQ DGVDVLAKVG GFDLVGMAGV MLGAASCGLP VLLDGFLSYA AALAACQMSP AIKPYLIPSH LSAEKGARIA LSHLGLEPYL NMEMRLGEGS GAALAMPIIE AACAIYNNMG ELAASNIVLP GNTTSDLNS
<b>Source</b>	Baculovirus
<b>Target Names</b>	cobT
<b>Protein Names</b>	Recommended name: Nicotinate-nucleotide--dimethylbenzimidazole phosphoribosyltransferase Short name= NN:DBI PRT EC= 2.4.2.21 Alternative name(s): N(1)-alpha-phosphoribosyltransferase
<b>Expression Region</b>	1-359
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.