



# Recombinant Escherichia coli O139:H28 Peptidase B (pepB)

<b>Product Code</b>	CSB-EP421964EJD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A7ZPW6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O139:H28 (strain E24377A / ETEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTEAMKITLS TQPADARWGE KATYSINNDG ITLHLNGADD LGLIQRAARK IDGLGIKHVQ LSGEGWDADR CWAFFWQGYKA PKGTRKVEWP DLDDAQRQEL DNRLMIIDWV RDTINAPAE LGPSQLAQRA VDLISNVAGD RVTYRITKGE DLREQGYMGL HTVGRGERS PVLLALDYNP TGDKEAPVYA CLVGKGITFD SGGYSIKQTA FMDSMKSDMG GAATVTGALA FAITRGLNKR VKLFLCCADN LISGNAFKLG DIITYRNGKK VEVMTDAEG RLVLADGLID ASAQKPELII DAATLTGAAK TALGNDYHAL FSFDDALAGR LLASASQENE PFWRLPLAEF HRSQLPSNFA ELNNTGSAAY PAGASTAAGF LSHFVENYQQ GWLHIDCSAT YRKAPVEQWS AGATGLGVRT IANLLTA
<b>Source</b>	E.coli
<b>Target Names</b>	pepB
<b>Protein Names</b>	Recommended name: Peptidase B EC= 3.4.11.23 Alternative name(s): Aminopeptidase B
<b>Expression Region</b>	1-427
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.