



Recombinant Escherichia coli O139:H28 UPF0597 protein yhaM (yhaM)

Product Code	CSB-EP421975EJD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZS00
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Purity	>85% (SDS-PAGE)
Sequence	MFDSTLNPLW QRYILAVQEE VKPALGCTEP ISLALAAAVA AAELEGPVER VEAWVSPNLM KNGLGVTVPG TGMVGLPIAA ALGALGGNAN AGLEVLKDAT AQAIADAKAL LAAGKVSVKI QEPCNEILFS RAKVWNGEKW ACVTIVGGHT NIVHIETHNG VVFTQQACVA EGEQESPLTV LSRTTLAEIL KFNVEVPFAA IRFILDSAKL NCALSQEGLS GKWGLHIGAT LEKQCERGLL AKDLSSSIVI RTSAASDARM GGATLPAMSN SGSGNQGITA TMPVVVVAEH FGADDERLAR ALMLSHLSAI YIHNQLPRLS ALCAATTAAM GAAAGMAWLV DGRYETISMA ISSMIGDVSG MICDGASNSC AMKVSTSASA AWKAVLMALD DTAVTGNEGI VAHDVEQSIA NLCALASHSM QQTDRQIIEI MASKAR
Source	E.coli
Target Names	yhaM
Protein Names	Recommended name: UPF0597 protein yhaM
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.