



Recombinant Escherichia coli O9:H4 Succinyl-diaminopimelate desuccinylase (dapE)

Product Code	CSB-MP422066EJF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8A2W3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	>85% (SDS-PAGE)
Sequence	MSCPVIELTQ QLIRRPSLSP DDAGCQALLI ERLQAIGFTV ERMDFADTQN FWAWRGQGET LAFAGHTDVV PPGDADRWIN PFEPTIRDG MLFGRGAADM KGSLAAMVVA AERFVAQHPN HTGRLAFLIT SDEEASAHNG TVKVVEALMA RNERLDYCLV GEPSSIEVVG DVVKNRRRS LTCNLTIHGV QGHVAYPHLA DNPVHRAAPF LNELVAIEWD QGNEFFPATS MQIANIQAGT GSNNVIPGEL FVQFNFRFST ELTDEMIKAQ VLALLEKHQL RYTVDDWWLSG QPFLTARGKL VDAVVNAVEH YNEIKPQLLT TGGTSDGRFI ARMGAQVVEL GPVNATIHKI NECVNAADLQ LLARMYQRIM EQLVA
Source	Mammalian cell
Target Names	dapE
Protein Names	Recommended name: Succinyl-diaminopimelate desuccinylase Short name= SDAP desuccinylase EC= 3.5.1.18 Alternative name(s): N-succinyl-L-lysine desuccinylase
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.