



Recombinant Escherichia coli O9:H4 Large ribosomal subunit protein uL4 (rpID)

Product Code	CSB-EP422085EJF-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8A5C4
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	>85% (SDS-PAGE)
Sequence	MELVLKDAQS ALTVSETTFG RDFNEALVHQ VVVAYAAGAR QGTRAQKTRA EVTGSGKKPW RQKGTGRARS GSIKSPIWRS GGVTFARPQ DHSQKVNKKM YRGALKSILS ELVRQDRLIV VEKFSVEAPK TKLLAQKLKD MALEDVLIIT GELDNLFLA ARNLHKVDVR DATGIDPVSL IAFDKVVMTA DAVKQVEEML A
Source	E.coli
Target Names	rpID
Protein Names	Recommended name: 50S ribosomal protein L4
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.