



# Recombinant Human Xaa-Pro dipeptidase (PEPD)

<b>Product Code</b>	CSB-EP017784HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P12955
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AAATGPSFW LGNETLKVPL ALFALNRQRL CERLRKNPAV QAGSIVVLQG GEETQRYCTD TGVLFRQESF FHWAFGVTEP GCYGVIVDVT GKSTLFPVRL PASHATWVGK IHSKEHFKEK YAVDDVQYVD EIASVLTSQK PSVLLTLRGV NTDSGSVCRE ASFDGISKFE VNNTILHPEI VECRVFKTDM ELEVLRNTK ISSEAHREVM KAVKVGMEKEY ELESLEFYHC YSRGGMRHSS YTCICGSGEN SAVLHYGHAG APNDRTIQNG DMCLFDMGGE YYCFASDITC SFPANGKFTA DQKAVYEAVL RSSRAVMGAM KPGVWWPDMH RLADRIHLEE LAHMGILSGS VDAMVQAHLG AVFMPHGLGH FLGIDVHDVG GYPEGVERID EPGLRSLRTA RHLQPGMVLV VEPGIYFIDH LLDEALADPA RASFLNREVL QRFRGFGGVR IEEDVVVTDS GIELLTCVPR TVEEIEACMA GCDKAFTPFSS GPK
<b>Source</b>	E.coli
<b>Target Names</b>	PEPD
<b>Protein Names</b>	Recommended name: Xaa-Pro dipeptidase Short name= X-Pro dipeptidase EC= 3.4.13.9 Alternative name(s): Imidodipeptidase Peptidase D Proline dipeptidase Short name= Prolidase
<b>Expression Region</b>	2-493
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.