



Recombinant Escherichia coli O9:H4 Peptidase T (pepT)

Product Code	CSB-EP423082EJF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZZ82
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	>85% (SDS-PAGE)
Sequence	MDKLLERFLN YVSLDTQSKA GVRQVPSTEG QWKLLHLLKE QLEEMGLINV TLSEKGTLM TLPANVPGDI PAIGFISHVD TSPDCSGKNV NPQIVENYRG GDIALGIGDE VLSPVMFPVL HQLLGQTLIT TDGKTLLGAD DKAGIAEIMT ALAVLQQKNI PHGDIRVAFT PDEEVGKGA K HFDVDAFDAR WAYTVDGGGV GELEFENFNA ASVNIKIVGN NVHPGTAKGV MVNALS LAAR IHA E VPADES PEMTEGYEGF YHLASMKGTV DRADMHYIIR DFDRKQFEAR KRKMMEIAKK VGKGLHPDCY IELVI EDSY Y NMREKVVEHP HILDIAQQAM RDCDIEPELK PIRGGTDGAQ LSFMGLPCPN LFTGGYNYHG KHEFVTLEGM EKAVQVIVRI AELTAQRK
Source	E.coli
Target Names	pepT
Protein Names	Recommended name: Peptidase T EC= 3.4.11.4 Alternative name(s): Aminotripeptidase Short name= Tripeptidase Tripeptide aminopeptidase
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.