



# Recombinant Escherichia coli O9:H4 Peptidase T (pepT)

<b>Product Code</b>	CSB-YP423082EJF
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A7ZZ82
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O9:H4 (strain HS)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDKLLERFLN YVSLDTQSKA GVRQVPSTEG QWKLLHLLKE QLEEMGLINV TLSEKGLTMA TLPANVPGDI PAIGFISHVD TSPDCSGKNV NPQIVENYRG GDIALGIGDE VLSPVMFPVL HQLLGQTLIT TDGKTLLGAD DKAGIAEIMT ALAVLQQKNI PHGDIRVAFT PDEEVGKGAK HFDVDAFDAR WAYTVDGGGV GELEFENFNA ASVNIKIVGN NVHPGTAKGV MVNALS LAAR IHAEPADES PEMTEGYEGF YHLASMKGTV DRADMHYIIR DFDRKQFEAR KRKMMEIAKK VGKGLHPDCY IELVIEDSY NMREKVVEHP HILDIAQQAM RDCDIEPELK PIRGGTDGAQ LSFMGLPCPN LFTGGYNYHG KHEFVTLEGM EKAVQVIVRI AELTAQRK
<b>Source</b>	Yeast
<b>Target Names</b>	pepT
<b>Protein Names</b>	Recommended name: Peptidase T EC= 3.4.11.4 Alternative name(s): Aminotripeptidase Short name= Tripeptidase Tripeptide aminopeptidase
<b>Expression Region</b>	1-408
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.