



# Recombinant Escherichia coli O139:H28 Ribulokinase (araB)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP424012EJD-B   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | A7ZHF4  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O139:H28 (strain E24377A / ETEC)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MAI AIGLDFG SDSVRALAVD CASGEEIATS VEWYPRWQKG QFCDAPNNQF<br>RHHPRDYIES MEAALKTVLA ELSVEQRAAV VGIGVDTTGS TPAPIDANGN<br>VLALRPEFAE NPNAMFVLWK DHTAVEEAE ITRLCHAPGN VDYSRYIGGI<br>YSSEWFWAKI LHVTRQDSAV AQSAASWIEL CDWVPALLSG TTRPQDIRRG<br>RCSAGHKSLW HESWGGLPPA SFFDELDPIL NRHLPSPLFT DTWTADIPVG<br>TLCPEWAQRL GLPESVVISG GAFDCHMGAV GAGAQPINALV KVI GTSTCDI<br>LIADKQSVGE RAVKGICGQV DGSVVPFGFIG LEAGQSAFGD IYAWFGRVLG<br>WPLEQLAAQH PELKAQINAS QKQLLPALTE AWAKNPSLDH LPVVDWVFN<br>R RTPNANQRL KGVITDLNLA TDAPLLFGGL IAATAFGARA IMECFTDQGI<br>AVNNVMALGG IARKNQVIMQ ACCDVLNRPL QIVASDQCCA LGAAIFAAVA<br>AKVHADIPSA QKMASAVEK TLQPRSEQAQ RFEQLYRRYQ<br>QWAMSAEQHY LPTSAPAQAA QAVATL |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | araB  |
| <b>Protein Names</b>     | Recommended name: Ribulokinase EC= 2.7.1.16   |
| <b>Expression Region</b> | 1-566   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.