



Recombinant Escherichia coli O139:H28 Methionine--tRNA ligase (metG)

Product Code	CSB-EP424049EJD-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZNT3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Purity	>85% (SDS-PAGE)
Sequence	<p>MTQVAKKILVTCALPYANGSIHLGHMLEHIQADVWVRYQRMRGHEVNFICADD AHGTPIM LKAQQLGITPEQMIGEMSQEHQTDFAGFNISYDNYHSTHSEENRQLSELIYSRL KENGF KNRTISQLYDPEKGMFLPDRFVKGTCPKCKSPDQYGDNCEVCGATYSPTHELIE PKSVVSG ATPVMRDSEHFFFDLPSFSEMLQAWTRSGALQEYVANKMQEWFESGLQQW DISRDAPYFG FEIPNAPGKYFYVWLDAPIGYMGSFKNLCDKRGDSVSFDEYWKKDSTAELYH FIGKDIVY FHSLFWPAMLEGSNFRKPTNLFVHGYVTVNGAKMSKSRGTFIKASTWLNHFD ADSLRYYY TAKLSSRIDDLNLEDFVQRVNADIVNKVVNLASRNAGFINKRFDGVLASELAD PQLYK TFTDAAEVIAGEAWESREFGKAVREIMALADLANRYVDEQAPWVVAKQEGRDA DLQAICSM GINLFRVLMTYLKPVLPLKTERAEAFNLTELTWDGIQQPLLGHKVNPFKALYNRI DMKQV EALVEASKEEVKAAAAPVTGPLADDPIQETITFDDFAKVDLRVALIENAEFVEGS DKLLR LTLDLGGEKRNVFSGIRSAYPDPQALIGRHTIMVANLAPRKMRFGISEGMVMA AGPGGKD IFLSPDAGAKPGHQVK</p>
Source	E.coli
Target Names	metG
Protein Names	Recommended name: Methionine--tRNA ligase EC= 6.1.1.10 Alternative name(s): Methionyl-tRNA synthetase Short name= MetRS
Expression Region	1-677
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.



Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.