



Recombinant Chicken Transforming protein p54/c-ets-1 (ETS1)

Product Code	CSB-EP007852CH-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P13474
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MKAAVDLKPT LTIKTEKVD IDLFPSPDME CADVPLLTPS SKEMMSQALK ATFSGFAKEQ QRLGIPKDPQ QWTETHVRDW VMWAVNEFSL KGVDFQKFCM NGAALCALGK ECFLELRPDF VGDILWEHLE ILQKEEAKPY PANGVNAAYP ESRYTSYFI SYGIEHAQCV PPSEFSEPSF ITESYQTLHP ISSEELLSLK YENDYPSVIL RDPVQTDSLQ TDYFTIKQEV VTPDNMCMGR ASRGKLGQD SFESIESYDS CDRLTQSWSS QSSFQSLQRV PSYDSFDSED YPAALPNHKP KGTFKDYVRD RADMNKDKPV IPAAALAGYT GSGPIQLWQF LLELLTDKSC QSFISWTGDG WEFKLSDPDE VARRWGKRKN KPKMNYEKLS RGLRYYYDKN IIHKTAGKRY VYRFVCDLQS LLGYTPEELH AMLDVKPDAD E
Source	E.coli
Target Names	ETS1
Protein Names	Recommended name: Transforming protein p54/c-ets-1
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.