



Recombinant Escherichia coli O9:H4 Flagellar L-ring protein (flgH)

Product Code	CSB-BP424124EJF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A7ZZ37
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	≥85% (SDS-PAGE)
Sequence	CAWIPSTPL VQGATSAQPV PGPTPVANGS IFQSAQPINY GYQPLFEDRR PRNIGDTLTI VLQENVASASK SSSANASRDG KTNFGFDTVP RYLQGLFGNA RADVEASGGN TFNGKGGANA SNTFSGLTV TVDQVLVNGN LHVVGKQIA INQGTEFIRF SGVVPRTIS GSNTVPSTQV ADARIEYVGN GYINEAQNMG WLQRFFLNLS PM
Source	Baculovirus
Target Names	flgH
Protein Names	Recommended name: Flagellar L-ring protein Alternative name(s): Basal body L-ring protein
Expression Region	22-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.