



Recombinant Escherichia coli O9:H4 N-acetylneuraminate epimerase (nanM)

Product Code	CSB-YP424192EJF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8A834
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	≥85% (SDS-PAGE)
Sequence	S VLPETPVPEK SGTGAIDNDT VYIGLGSGAGT AWYKLDTPQAK DKKWTALAAF PGGPRDQATS AFIDGNLYVF GGIGKNSKGL TQVFNDVHKY NPKTNSWVKL MSHAPMGMAH HVTFFVHNGKA YVTGGVNQNI FNGYFEDLNE AGKDSATIDK INAHYFDKKA EDYFFNKFLS SFDPTSTQQWS YAGESPWYGT AGAAVVKNGD KTWLVNGEAK PGLRTDAVFE LDFTGNLKW NKLAPVSSPD GVAGGFAGIS NDSLIFAGGA GFKGSRENYQ NGKNYAHEGL KKSYSTDIHL WHNGKWDKSG ELSQGRAYGV SLPWNNSLLI IGGETAGGKA VTDSVLISVK DNKVTVQN
Source	Yeast
Target Names	nanM
Protein Names	Recommended name: N-acetylneuraminate epimerase EC= 5.1.3.24 Alternative name(s): N-acetylneuraminate mutarotase Short name= Neu5Ac mutarotase Sialic acid epimerase
Expression Region	20-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.