



Recombinant Escherichia coli O139:H28 Enolase (eno)

| | |
|--------------------------|---|
| Product Code | CSB-EP425108EJD-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | A7ZQM2 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O139:H28 (strain E24377A / ETEC) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSKIVKIIGR EIIDSRGNPT VEAEVHLEGG FVGMAAAPSG ASTGSREALE LRDGDKSRFL GKGVTKAVAA VNGPIAQALI GKDAKDQAGI DKIMIDLDTG ENKSKFGANA ILAVSLANAK AAAAAKGMPL YEHIAELNGT PGKYSMPVPM MNIINGGEHA DNNVDIQEFM IQPVGAKTVK EAIRMGSEVF HHLAKVLKAK GMNTAVGDEG GYAPNLGSNA EALAVIAEAV KAAGYELGKD ITLAMDC AAS EFYKDGKYVL AGEGNKAFTS EEFTHFLEEL TKQYPIVSIE DGLDESDWDG FAYQTKVLGD KIQLVGDDL FVTNTKILKEG IEKGIANSIL IKFNQIGSLT ETLAAIKMAK DAGYTAVISH RSGETEDATI ADLAVGTAAG QIKTGSMRSRS DRVAKYNQLI RIEEALGEKA PYNGRKEIKG QA |
| Source | E.coli |
| Target Names | eno |
| Protein Names | Recommended name: Enolase EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 2-phosphoglycerate dehydratase |
| Expression Region | 1-432 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |