



# Recombinant Escherichia coli O9:H4 NADPH-dependent 7-cyano-7-deazaguanine reductase (queF)

<b>Product Code</b>	CSB-MP425194EJF
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A8A3S9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O9:H4 (strain HS)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSSYANHQAL AGLTLGKSTD YRDTYDASLL QGVPRSLNRD PLGLKADNLP FHGTDIWTLY ELSWLNAGKL PQVAVGHVEL DYTSVNLIES KSKFLYLNSF NQTRFNNWDE VRQTLERDLS TCAQGGKISVA LYRLDELEGQ PIGHFNGTCI DDQDITIDNY EFTTDYLENA TCGEKVVVEET LVSHLLKSNC LITHQPDWGS IQIQYRGRQI DREKLLRYLV SFRHHNEFHE QCVERIFNDL LRFCQPEKLS VYARYTRRGG LDINPWRSNS DFVPSTTRLV RQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	queF
<b>Protein Names</b>	Recommended name: NADPH-dependent 7-cyano-7-deazaguanine reductase EC= 1.7.1.13 Alternative name(s): 7-cyano-7-carbaguanine reductase NADPH-dependent nitrile oxidoreductase PreQ(0) reductase
<b>Expression Region</b>	1-282
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.