



Recombinant Escherichia coli O9:H4 1,4-alpha-glucan branching enzyme GlgB (glgB)

Product Code	CSB-EP425209EJF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8A5P2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Purity	>85% (SDS-PAGE)
Sequence	MSDRIDRDVI NALIAGHFAD PFSVLGMHKT TAGLEVRALL PDATDVWVIE PKTGRKLAKL ECLDSRGFFS GVIPRRKNFF RYQLAVVWHG QQNLIDDPYR FGPLIQEMDA WLLSEGTHLR PYETLGAHAD TMDGVTGTRF SVWAPNARRV SVVGQFNWD GRRHPMRLRK ESGIWELFIP GAHNGQLYKY EMIDANGNLR LKSDPYAFEAF QMRPETASLI CGLPEKVVT EERKKANQFD APISIEVHL GSWRRHTDNN FWLSYRELAD QLPYAKWVG FTHLELLPIN EHPFDGSWGY QPTGLYAPTR RFGTRDDFRY FIDAAHAAGL NVILDWVPGH FPTDDFALAE FDGTNLYEHS DPREGYHQDW NTLIYNYGRR EVSNFLVGNA LYWIERFGID ALRVDVASM IYRDYSRKEG EWIPNEFGGR ENLEAIEFLR NTNRILGEQV SGAVTMAEES TDFPGVSRPQ DMGGLGFWYK WNLGWMHDTL DYMKLDPVYR QYHHDKLTFG ILYNYTENFV LPLSHDEVVH GKKSILDRMP GDAWQKFANL RAYYGWMMWAF PGKLLFMGN EFAQGREWNH DASLDWHLLE GGDNWHHGVSQ RLVLDLNLTY RHHKAMHELD FDPYGFELV VDDKERSVLI FVRRDKEGNE IIVASNFTPV PRHDYRFGIN QPGKWREILN TDSMHHYHGSN AGNGGTVHSD EIASHGRQHS LSLTLPLPLAT IWLVREAE
Source	E.coli
Target Names	glgB
Protein Names	Recommended name: 1,4-alpha-glucan branching enzyme GlgB EC= 2.4.1.18 Alternative name(s): 1,4-alpha-D-glucan:1,4-alpha-D-glucan 6-glucosyl- transferase Alpha-(1->4)-glucan branching enzyme Glycogen branching enzyme Short name=
Expression Region	1-728
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.