



# Recombinant Escherichia coli O9:H4 3-octaprenyl-4-hydroxybenzoate carboxy-lyase (ubiD)

<b>Product Code</b>	CSB-EP425218EJF-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A8A6U8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O9:H4 (strain HS)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKYNDLRDFL TLLEQQGELK RITLPVDPHL EITEIADRTL RAGGPALLFE NPKGYSMPVL CNLFGTPKRV AMGMGQEDVS ALREVGKLLA FLKEPEPPKG FRDLFDKLPQ FKQVLNMPTK RLRGAPCQK IVSGDDVDLN RIPIMTCWPE DAAPLITWGL TVTRGPHKER QNLGIYRQQL IGKNKLIMRW LSHRGGALDY QEWCAAHPGE RFPVSVALGA DPATILGAVT PVPDTLSEYA FAGLLRGTKT EYVKCISNDL EVPASAEIVL EGYIDPGEMA PEGPYGDHTG YYNEVDSFPV FTVTHITQRE DAIYHSTYTG RPPDEPAVLG VALNEVFVPI LQKQFPEIVD FYLPPEGCSY RLAVVTIKKQ YAGHAKRMMM GVWSFLRQFM YTKFVIVCDD DVNARDWWDV IWAITTRMDP ARDTVLVENT PIDYLDFA SP VSGLGSKMGL DATNKWPGET QREWGRPIKK DPDVVAHIDA IWDELAIFNN GKSA
<b>Source</b>	E.coli
<b>Target Names</b>	ubiD
<b>Protein Names</b>	Recommended name: 3-octaprenyl-4-hydroxybenzoate carboxy-lyase EC= 4.1.1.- Alternative name(s): Polyprenyl p-hydroxybenzoate decarboxylase
<b>Expression Region</b>	1-494
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.