



Recombinant Rat NADP-dependent malic enzyme (Me1)

Product Code	CSB-EP013632RA-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P13697
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MDPRAPRRRH THQRGYLLTR DPHLNKDLAF TLEERQQLKI HGLLPPCIVN QEIQVLRVIK NFERLNSDFD RYLLMLDLQD RNEKLFYSVL MSNVEKFMPI VYTPTVGLAC QQYSLAFRKP RGLFISIHDK GHIASVLNAW PEDVVKAIVV TDGERILGLG DLGCNGMGIP VGKLALYTAC GGVNPQQCLP ITLDVGTENE ELLKDPLYIG LRHRRVRGPE YDAFLDEFME AASSKYGMNC LIQFEDFANL NAFRLLNKYR NKYCTFNDDI QGTASVAVAG LLAALRITKN KLSDQTVLFQ GAGEAALGIA HLIVMAMEKE GLSKEKARQK IWLVDKGLI VKGRASLTEE KEVFAHEHEE MKNLEAIVQK IKPTALIGVA AIGGAFTEQI LKDMAAFNER PIIFALSNT SKAECSAEEC YKVTKGRAIF ASGSPFDPVT LPDGRTLFP QGNNSYVFPG VALGVVACGL RHINDSVFLT TAEVISQQVS DKHLEEGRLY PPLNTIRDVS LKIAVKIVQD AYKEKMATVY PEPQNKEEFV SSQMYSTNYD QILPDCYSWP EEVQKIQT KV NQ
Source	E.coli
Target Names	Me1
Protein Names	Recommended name: NADP-dependent malic enzyme Short name= NADP-ME EC= 1.1.1.40 Alternative name(s): Malic enzyme 1
Expression Region	1-572
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes a cytosolic, NADP-dependent enzyme that generates NADPH for fatty acid biosynthesis. The activity of this enzyme, the reversible oxidative decarboxylation of malate, links the glycolytic and citric acid cycles. The regulation of expression for this gene is complex. Increased expression can result from elevated levels of thyroid hormones or by higher proportions of carbohydrates in the diet.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.