



# Recombinant Escherichia coli O139:H28 2- (5''-triphosphoribosyl)-3'-dephosphocoenzyme-A synthase

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP426062EJD   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | A7ZIZ9  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O139:H28 (strain E24377A / ETEC)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSMPATSTKT TKLATSLIDE YALLGWRAML TEVNLSPKPG LVDRINCGAH<br>KDMALED FHR SALAIQGWLP RFIEFGACSA EMAPEAVLHG LRPIGMACEG<br>DMFRATAGVN THKGSIFSLG LLCAAIGRLL QLNQPVTPTT VCSTAASFCR<br>GLTDRELRTN NSQLTAGQRL YQQLGLTGAR GEAEAGYPLV LNHALPHYLT<br>LLDQGLDPEL ALLDTLLLLM AINGDTNVAS RGEGGGLRWL QREAQTL LQK<br>GGIRTPADLD YLRQFDRECI ERNLSPGGSA DLLILT WFLA QI   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | citG  |
| <b>Protein Names</b>     | Recommended name: 2-(5''-triphosphoribosyl)-3'-dephosphocoenzyme-A synthase Short name= 2-(5''-triphosphoribosyl)-3'-dephospho-CoA synthase EC= 2.7.8.25  |
| <b>Expression Region</b> | 1-292   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |