



# Recombinant Escherichia coli O139:H28 DNA-directed RNA polymerase subunit alpha (rpoA)

<b>Product Code</b>	CSB-EP426112EJD-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A7ZSI4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O139:H28 (strain E24377A / ETEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQGSVTEFLK PRLVDIEQVS STHAKVTLEP LERGFHTLG NALRRILLSS MPGCAVTEVE IDGVLHEYST KEGVQEDILE ILLNLKGLAV RVQKGDEVIL TLNKSGIGPV TAADITHDGD VEIVKPQHVI CHLTDENASI SMRIKVQRGR GYVPASTRIH SEEDERPIGR LLVDACYSPV ERIAYNVEAA RVEQRTDLDK LVIEMETNGT IDPEEAIRRA ATILAEQLEA FVDLRDVRQP EVKEEKPEFD PILLRPVDDL ELTVRSANCL KAEAIHYIGD LVQRTEVELL KTPNLGKKSL TEIKDVLASR GLSLGMRLLEN WPPASIADE
<b>Source</b>	E.coli
<b>Target Names</b>	rpoA
<b>Protein Names</b>	Recommended name: DNA-directed RNA polymerase subunit alpha Short name= RNAP subunit alpha EC= 2.7.7.6 Alternative name(s): RNA polymerase subunit alpha Transcriptase subunit alpha
<b>Expression Region</b>	1-329
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.