



Recombinant Brugia malayi tRNA (guanine-N (7)-)-methyltransferase

Product Code	CSB-MP427056BWV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8NFF0
Product Type	Recombinant Protein
Immunogen Species	Brugia malayi (Filarial nematode worm)
Purity	>85% (SDS-PAGE)
Sequence	MVSTENKIGL FKNKDDDDIDG EEMRELPQKK FYRQRAHANP ISDHEFDYPV FPEQMDWKKY FGDFSEGRQV EFADVCGCYG GLLIKLSTLY PEALMVGLEI RVKVS DYVQD KIHALLRLREP GNYRNVACLR TNAMKYLPNY FRRHQLTKMF FLYPDPHFKK AKHKWRIITP TLLAEYAYVL KPGGLVYTIT DVEELHIWMV RHLSAHPLFE RLTDLEMKMD PVVEMLYDST EEGQKVARNE GSKWSAVFRR LPNPVLSS
Source	Mammalian cell
Target Names	Bm1_01445
Protein Names	Recommended name: tRNA (guanine-N(7)-)-methyltransferase EC= 2.1.1.33 Alternative name(s): tRNA (guanine(46)-N(7))-methyltransferase tRNA(m7G46)-methyltransferase
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.